

Bed Leaving Alarm with Voice Reminder and Pager

Code: MP-CTM-3BKIT

Price: £112.80 incl. VAT

Date: 01/12/2021

VAT Relief Eligible: Yes (£94.00 with VAT relief)

Delivery: £3.95 incl. VAT per order



- Part of the MPPL system
- Allows a carer to be alerted/alarmed when a person leaves their bed and can help reduce the risk of falls and accidents
- Kit includes everything required to get started, including batteries
- Record an audio message of up to 15 seconds, that's played every time someone leaves the bed (optional)

Bed Leaving Sensor Mat

- Thin and lightweight Bed Sensor detects exit from bed and sends alarm to pager
- 15-minute delayed alarm (optional) to accommodate bathroom visits etc
- Pager low battery alert signalled by red flashing LED, repeated every 2 weeks before battery failure
- Dimensions: 750 x 450mm

CTM3 Voice Alert Audio Alarm

- Supplied with slot-in stand for easy mounting
- Recording triggered when movement from the bed is detected
- Repeats for 5 minutes, or until dismissed at the alarm
- If pressure is re-applied, voice message halts and monitor resets
- Clear reset button for one-press recording cancellation
- Voice alert, tone only, or silent options
- Anti-tamper panel screws in place
- Low battery signalled by neon blue flash every 3 hours, repeated for 2 weeks before battery failure
- Nurse call relay output option for care homes
- Supplied with 3 x AAA batteries

What if I don't want the voice alert function?

The voice alert audio alarm serves two purposes;

1. Transmits the alarm from the pressure mat to the pager
2. Plays a prerecorded audio message

Due to no.1, the voice alert audio alarm is a required component of this kit, however, if you do not require an audio message to be played there is a silent option available for the voice alert function.

MPPL Pager

- Supplied with 2 x AA batteries
- Detailed information can be found on our [MPPL pager page](#)

Contents

- 1 x Bed leaving sensor mat
- 1 x A voice alert audio alarm
- 1 x MPPL pager

Additional Images



How to Order

Orders can be placed securely online at AlzProducts.co.uk - or over the phone with our UK sales team.

Call 024 7642 2224 to place your order or for further information. This product's order code is MP-CTM-3BKIT.