

Chair Leaving Alarm with LCD Pager

Code: 1005415

Price: £121.20 incl. VAT

Date: 15/06/2019

VAT Relief Eligible: No

Delivery: £3.95 incl. VAT per order



Our **Chair Leaving Alarm with LCD Pager** system provides alerts for chair rise with an alarm transmitted to the new Pager with LCD display.

The pager will link to all Medpage transmitter products. Ideal for home use and professional care. The Pager can also be used with door alarms, PIR, and floor mats.

Sensors are manufactured in high quality polyvinyl and are fitted with a 2 Metre lead for connection to the monitor transmitter. Normal strength disinfectants can be used to sterilise the sensors if required.

Chair sensor Measures 15" x 12"

LCD Pager:

- Pager featured with vibration only or tone beeper alert
- Pager has LCD display for device ID, location and general system status
- Operates with any Medpage or Easylink transmitter products
- Powered by 2 x AA batteries
- Pager Dimensions: W 55mm H 75mm D 25mm
- Strong belt clip
- Low battery warning
- Alarm volume control
- Alarm display, volume control, stored alarm review button
- Radio integrity monitoring
- Multi-User alarm identification display allows up to 99 devices to signal pager

Chair Sensor Transmitter:

- 15 minute alarm delay option
- Sensor Dimensions: H 115mm W 75mm D 20mm
- Case material ABS
- Wall fixing lugs on rear of case
- Transmission frequency 433.92 MHz
- Range 200M (open field)
- Power requirement 9V PP3 or mains 12V DC power adaptor
- Low battery warning sound + flashing LED
- Built in nurse call button, alarm reset button

How to Order

Orders can be placed securely online at AlzProducts.co.uk - or over the phone with our UK sales team.

Call 024 7642 2224 to place your order or for further information. This product's order code is 1005415.