

Chair Occupancy Alarm Mat with Transmitter

Code: 1005464

Price: £76.80 incl. VAT

Date: 28/07/2020

VAT Relief Eligible: Yes (£64.00 with VAT relief)

Delivery: £3.95 incl. VAT per order



The **Chair Occupancy Alarm Mat with Transmitter** is an ideal way of monitoring a person's movements. Positioned on the chair, the sensor sounds an alert when weight is removed.

- Part of the MPPL system
- Allows a carer to be alerted/alarmed when a person leaves their chair and can help reduce the risk of falls and accidents
- Record an audio message of up to 15 seconds, that's played every time someone leaves the chair (optional)

Chair Leaving Sensor Mat

- Thin and lightweight Chair Sensor that detects when someone leaves the chair and sends an alarm to the pager
- 15-minute delayed alarm (optional) to accommodate bathroom visits etc
- Pager low battery alert signalled by red flashing LED, repeated every 2 weeks before battery failure
- Dimensions of 380 x 300mm

CTM3 Voice Alert Audio Reminder

- Supplied with slot-in stand for easy mounting
- Recording triggered when movement from the chair is detected
- Repeats for 5 minutes before going into power saving mode
- If pressure is re-applied, voice message halts and monitor resets
- Clear reset button for one-press recording cancellation
- Voice alert, tone only, or silent options
- Anti-tamper panel screws in place
- Low battery signalled by neon blue flash every 3 hours, repeated for 2 weeks before battery failure
- Nurse call relay output option for care homes

What if I don't want the voice alert function?

The voice alert audio alarm serves two purposes;

1. Transmits the alarm from the pressure mat to the pager
2. Plays a prerecorded audio message

Due to no.1, the voice alert audio alarm is a required component of this kit, however, if you do not require an audio message to be played there is a silent option available for the voice alert function.

Contents

- 1 x Chair leaving sensor mat
- 1 x A voice alert audio alarm

How to Order

Orders can be placed securely online at AlzProducts.co.uk - or over the phone with our UK sales team.

Call 024 7642 2224 to place your order or for further information. This product's order code is 1005464.