

Bed/Chair Exit Alarm with Magnetic Clothing Clip Removal Alarm


Code: 1005471

Price: £74.50 incl. VAT

Date: 08/11/2019

VAT Relief Eligible: No

Delivery: £3.95 incl. VAT per order



This **Bed/Chair Exit Alarm with Magnetic Clothing Clip Removal Alarm** offers state-of-the-art fall prevention technology for home and professional care use.

The dual alarm trigger function makes it really quick and easy to set up.

Clothing Clip - secure the sensor alarm controller to the bed head, armchair or wheelchair. Fasten the clip to the persons clothing, attach the magnet to the alarm sensor, adjust the length of the cord and its ready to go.

When the person moves and pulls the magnet away from the controller, the selected alarm option activates.

Features

- Personal message recording up to 15 seconds to play during alarm "Please stay in bed"
- Tone alarm with volume control (Switch selectable)
- Microphone jack input to upload MP3 to play during alarm
- Reset button on monitor face
- Pull string magnet alarm (attach to clothing)
- Port to connect with nurse call system (NC/NO) 1/4" jack socket (6.3mm)
- Operates on batteries or via AC adaptor
- Low battery and nurse call plug removal alert
- Pad lost alert

Technical Specification

- Power Supply: 6V AC adaptor, 300mA or 3 x AA-type alkaline batteries
- Battery life: 1 year (Nominal)
- Dimensions: (L x W x D): 120 x 77 x 34mm
- Housing material: ABS (UL-94HB, flammability grade) Net Weight: 135g (without soft cover & batteries)
- Operating temp. 5°C ~ 40°C
- Humidity tolerance < 85%

Included

- 1 x BMDV-11 sensor alarm
- 1 x Protective rubber cover
- 1 x Clothing clip/alarm magnet
- 3 x AA batteries
- 1 x Protective rubber sleeve

How to Order

Orders can be placed securely online at AlzProducts.co.uk - or over the phone with our UK sales team.

Call 024 7642 2224 to place your order or for further information. This product's order code is 1005471.