

Chair Leaving Alarm with Nurse Call Connection

Code: 1005583

Price: £74.40 incl. VAT

Date: 06/07/2020

VAT Relief Eligible: Yes (£62.00 with VAT relief)

Delivery: £3.95 incl. VAT per order



This new advanced **Chair Leaving Alarm with Nurse Call Connection** system provides alerts when someone leaves their chair.

Sensors are manufactured in high quality polyvinyl and are fitted with a 2 Metre lead for connection to the monitor transmitter. Normal strength disinfectants can be used to sterilise the sensors if required.

Chair Sensor Measures: 120 x 77 x 34mm

Chair Sensor Transmitter

- Case material ABS
- Wall fixing lugs on rear of case
- Transmission frequency 433.92 MHz
- Range 200M (open field)
- Power requirement 9V PP3 or mains 12V DC power adaptor
- Low battery warning sound + flashing LED
- Built in nurse call button, alarm reset button

How to Order

Orders can be placed securely online at AlzProducts.co.uk - or over the phone with our UK sales team.

Call 024 7642 2224 to place your order or for further information. This product's order code is 1005583.